

# Infection control and patient care.

They're one  
and the  
same.

# Witness efficient *Clostridioides difficile* detection in as little as 24 hours!



Hardy Diagnostics introduces C diff Banana Broth™, which is recommended as a selective detection broth for the recovery of *Clostridioides difficile* vegetative cells or spores from environmental samples. 20/pk **K226**

- Specificity found to be 100%. Ask for reference
- No special anaerobic incubation equipment required
- Selective and differential medium for recovering *Clostridioides difficile* spores and vegetative cells
- No need to pre-reduce the media
- Easy color change read-out



Negative

Positive

## Verify the efficacy of your cleaning protocol

# Accessory Products

## Saline

For pre-moistening of the swab.  
0.85%, 2ml, 13x100mm tube, shatter  
resistant, polycarbonate tube.  
20/pkg

R45



## HydraFlock

Flocked nylon swab  
tips act like a soft brush and allow  
better collection and elution samples.  
Swabs are scored, sterile, and  
individually wrapped.  
10 boxes of 50

FS1HD



## KWIK-STIK™ (Positive Control)

A ready-to-use swab that contains a single  
lyophilized pellet of a single microorganism strain,  
*Clostridioides difficile*, ATCC® 9689™. To use, simply  
crush the rehydrating ampule, mix, remove swab,  
and inoculate, by Microbiologics®.  
6/pkg

0329K







## Boekel Incubator

- Overall dimensions:  
13"x14"x16-1/4"
- Temp. range:  
ambient +5° C to  
60° C
- LCD temperature  
indicator
- Supplied with two  
removable zinc  
plated steel shelves  
Each 133000



# Environmental detection of *Clostridioides difficile*



Evaluation samples  
available!



**Made in the USA**

### Headquarters

1430 West McCoy Lane  
Santa Maria, CA 93455  
800.266.2222 : phone  
Sales@HardyDiagnostics.com  
HardyDiagnostics.com

### Distribution Centers

Santa Maria, California  
Olympia, Washington  
Salt Lake City, Utah  
Phoenix, Arizona  
Dallas, Texas  
Springboro, Ohio  
Lake City, Florida  
Albany, New York  
Raleigh, North Carolina

The Quality Management System at the Hardy Diagnostics manufacturing facility is certified to ISO 13485.



ISO 13485  
FM 572526

