




Is *C diff* lurking
in your hospital?

 Online Ordering Available

C DIFF BANANA BROTH

Witness Efficient Clostridium difficile Detection!
In as little as 24 hours!



Verify the efficacy of your cleaning protocol.

C diff Banana Broth™

Hardy Diagnostics introduces C diff Banana Broth™, which is recommended as a selective detection broth for the recovery of *Clostridium difficile* vegetative cells or spores from environmental samples.

20/pk

K226



ACCESSORY PRODUCTS

Saline

For pre-moistening of the swab. 0.85%, 2ml, 13x100mm tube, shatter resistant, polycarbonate tube.

20/pk

R45



HydraFlock



Flocked nylon swab tips act like a soft brush and allow better collection and elution samples. Swabs are scored, sterile, and individually wrapped, by Hardy Diagnostics.

10 boxes of 50

FS1HD

Incubator


- Overall dimensions: 13"x14"x16-1/4"
- Temp. range: ambient +5° C to 60° C
- LCD temperature indicator
- Supplied with two removable zinc plated steel shelves



Small Counter-top Incubator

133000

ENVIRONMENTAL DETECTION OF C. DIFFICILE



Evaluation samples available!

KWIK-STIK™ (Positive Control)

A ready-to-use swab that contains a single lyophilized pellet of a single microorganism strain. *Clostridium difficile*, ATCC® 9689™. To use, simply crush the rehydrating ampule, mix, remove swab, and inoculate, by Microbiologics®.

6/pk

0329K



A Culture of Service™

1430 West McCoy Lane

Santa Maria, CA 93455

www.HardyDiagnostics.com | Sales@HardyDiagnostics.com

Distribution Centers: California | Washington | Utah | Arizona

Ohio | Texas | New York | Florida | North Carolina