



HARDY DIAGNOSTICS

Dilu-Lok II™

The Dilu-Lok II™ eliminates the time and money spent on preparing your own dilution vials! Our pre-filled dilution vial is designed especially for the food, dairy, pharmaceutical, and water industries. The vials are filled to 25ml, 90ml, or 99ml with a variety of buffering and enrichment solutions. Each lot is meticulously tested for pH, sterility, fill volume, and toxicity.

NEW! Call to receive a FREE sample pack to evaluate in your lab!

Save time and money with the Dilu-Lok II™ pre-filled dilution vial



Simple...
Easy-to-remove, perforated, shrink wrap ensures product integrity.



One-Handed...
Hinged lid allows for easy opening.



Large Opening...
A wide 45mm opening accommodates most sample sizes.



Leak-Proof...
Recently re-engineered for an improved seal.

Please order in multiples of 50.

Product	Cat. no.
Butterfield's Phosphate, 99ml, 50 per box	D599
Butterfield's Phosphate, 90ml, 50 per box	D590
Butterfield's Phosphate, 25ml, 50 per box	D125
Peptone Water 0.1%, 99ml, 50 per box	D299
Peptone Water 0.1%, 90ml, 50 per box	D290
Phosphate Buffer with MgCl, 99ml, 50 per box ...	D699
Phosphate Buffer with MgCl, 90ml, 50 per box ...	D690



Santa Maria, CA | Salt Lake City, UT | Phoenix, AZ | Tacoma, WA
 Lima, OH | Dallas, TX
 Phone: 800.266.2222
 Fax: 805.346.2760
 Website: www.HardyDiagnostics.com
 Email: Sales@HardyDiagnostics.com