

Products for Influenza

 **HARDY**
DIAGNOSTICS
A Culture of Service™

 Online Ordering Available

Chlamydia, Viral, Mycoplasma Collection & Transport

TransPRO[™] CVM Transport System is a specimen collection kit featuring PurFlock[®] Ultra swab and room temperature stability.

- Ready-to-use transport system
- Premium medical grade plastic components to improve product shelf life
- Advanced kit design for reliability and ease-of-use
- Stable medium at room temperature
- Soft tip for patient comfort
- Economically priced
- Flexible shaft

Ordering Information

Specimen Collection Kit with PurFlock[®] Ultra Swab

Adult/Pediatric,
3ml, 50/pk.....**TPCVM1**

TransPRO[™] CVM Transport System is a specimen collection kit featuring the PurFlock Ultra[®] swab and transport system. It is the ideal choice with excellent absorption and elution of specimens. For collection and preservation of Virus, Chlamydia, Mycoplasma, and Ureaplasma.

The PurFlock Ultra[®] swab is the optimal swab for various rapid diagnostic tests including specimen collection, buccal cell collection, DNA testing, forensic evidence collection, and more.

A Comparison of Features

Main Competitive Advantages in Favor of TransPRO[™] CVM

TransPRO [™] CVM by Hardy Diagnostics	UTM [™]	
Pack Size:	50/pk	300/pk
Intended Use	Viruses, chlamydia, mycoplasma, ureaplasma	Viruses, chlamydia, mycoplasma, ureaplasma
Storage Temperature	2-25 °C	2-25 °C
Transport Temperature	4-25 °C for 48 hours	4-25°C for 48 hours
Swab Type	Naso-pharyngeal	Naso-pharyngeal
Swab-Cap Capture	Yes	Standard tip only
Swab Tip	PurFlock [®] Ultra (Patented New Technology)	Nylon Flocked (Patented)
Swab Tip Color	White	Yellow tinge
Cap Color	Red	Red
Contains Gelatin	Yes	Yes
Manufacture Location	USA	Italy
Shelf Life	18 months	12 months
Fill Volume	3ml	3ml
Cat. no.	TPCVM1	346C

Swab-Cap Capture Feature:

After the sample is collected, the swab is placed into the CVM Medium and the shaft is snapped at the score point. The lid is then screwed on. When the lid is unscrewed, the swab is connected to the lid.

Polystyrene Shaft

Compatible with most automated systems.



PurFlock® Ultra Swab Tip Technology for maximum specimen recovery.
Clean, white, effective!

Vial: Overall Length
16 x 100 mm

TransPRO™ CVM Medium



Score Point

Call 800.266.2222 for a free evaluation sample!

Adult/Pediatric, 3ml, 50/pk.....**TPCVM1**

Packaging includes:
Vial with 3ml CVM Medium
and one PurFlock® Ultra
swab, individually wrapped.



Made in the USA

Collection Swabs

Flocked Swab (nylon), Flexible Shaft, Mini-tip, Nasopharyngeal (NP)

This nylon-flocked swab is a flexible, thin, plastic shaft for Nasopharyngeal specimen collection, scored at 123mm from the tip, individually wrapped, and sterile.

500 swabs/caseFS3HD

Dacron Swab, Nasopharyngeal (NP)

Sterile, straight aluminum wire; dacron tipped swab, single, and individually wrapped.

100 per box.....4175BX


1000 per case4175CS

Rayon Swab, Transporter Nasopharyngeal (NP) with tube

A sterile rayon swab is on a twisted wire shaft, very flexible, in a plastic transport tube, without transport media. Individually wrapped.

100/box4168BX

1000/box4168CS



Mini-Tip, Nylon Flocked Swab

Influenza Test Kits





QuickVue™ Influenza A+B kit

Tests for and differentiates between Influenza A or B, 30 seconds of hands-on time, 10 minutes total time, and easy read-out. CLIA waived.

25 tests kit.....20183

QuickVue™ Influenza A+B Kit



Hand Sanitizers

Smart-San™ - Hand Sanitizer Spray

This ethanol-based spray with emollients keeps hands soft. Fragrance free, rated to kill 99.99% of Gram-negative and positive bacteria. Requires pump, sold separately.

Six/one liter bottles SMA0016U
 Pump JPP10007



Smart-San™ Hand Sanitizer Spray



AloeSafe™ - Hand Sanitizer Gel

This alcohol gel with aloe vera, vitamin E, and lavender will not dry hands, keeps skin soft, and conforms to CDC guidelines for hand sanitizers.

- Contains healing Vitamin E and moisturizing *Aloe vera*
- Does not cause drying
- Unique formula that's effective without harsh ingredients
- Does not contain propylene glycol (automobile antifreeze)
- Does not contain sticky glycerin
- Does not contain estrogenic parabens
- Kills 99.9% of bacteria

2oz pocket sizeAS2HD
 4oz travel sizeAS4HD
 16oz pump bottleAS16HD
 Wall mounted 16oz bottle holderASH16



AloeSafe™ Hand Sanitizer Gel

Face Masks

Masks, N95, Particulate Respirator

The N95 particulate respirator mask meets WHO and CDC recommended protection levels. It also meets NIOSH-defined bacterial filtration efficacy against airborne pathogens, and contamination.

- ✓ Mask held away from the nose and mouth for comfortable breathing with Magic Arch® feature
- ✓ Positive Facial Lock® with the Twist-Seal® chin piece design helps for a perfect fit with no gaps that might allow small particles to be inhaled around the edges of the mask
- ✓ Downward vision not restricted for added safety
- ✓ NIOSH-approved masks
- ✓ One size fits all

210 masks/case..... USP10206



N95, Particulate Respirator Mask

Disinfectants


CaviCide™ Disinfectant

This best all-around disinfectant contains quarternary ammonium and 17% alcohol with detergent and is recommended for benchtop cleaning and disinfection. Non-toxic and ready-to-use, it has a broad spectrum of kill against bacteria, viruses, and fungi. It contains no harmful phenols, glutaraldehyde, or bleach.

2oz travel size.....	131002
8oz spray bottle.....	131008
24oz bottle with sprayer	131024
One gallon.....	131000



*CaviCide™ Disinfectant,
Cat. no. 131000.*





Clorox® Broad Spectrum Quaternary, Disinfectant Cleaner

Fast-Acting Cleaner Disinfectant

- 33% faster than a leading Quaternary Disinfectant*
- Only quaternary disinfectant registered to kill representatives of all ESKAPE pathogens*, which cause two-thirds of all HAIs
- Kills MRSA, HIV, TB, VRE, Hepatitis B, C, and more
- Kills bacteria in 2 minutes
- Fragrance-free
- Gentle on surfaces

Healthlink, 32oz30649

*Based on Federal Master Label



Pop-up Disinfectant Wipes

DisCide™ Pop-up Wipes

These disinfecting towelettes are impregnated with a second-generation quaternary ammonium and isopropyl alcohol. The disinfectant towelette kills TB, HIV (AIDS) virus, Herpes 2, and Adenovirus 2 in one minute. Towelettes offer user-friendly, non-aerosol application in the laboratory.

60 per container.....10DIS
 160 per container.....60DIS
 Wall rack for above.....109

HARDY DIAGNOSTICS
 A Culture of Service™

1430 WEST MCCOY LANE
 SANTA MARIA, CA 93455

CUSTOMER SERVICE
 PHONE: 800.266.2222
 FAX: 805.346.2760
 SALES@HARDYDIAGNOSTICS.COM
 WWW.HARDYDIAGNOSTICS.COM

DISTRIBUTION CENTERS
 SANTA MARIA, CALIFORNIA
 LACEY, WASHINGTON
 SALT LAKE CITY, UTAH
 PHOENIX, ARIZONA
 DALLAS, TEXAS
 SPRINGBORO, OHIO
 LAKE CITY, FLORIDA
 ALBANY, NEW YORK
 NORTH CAROLINA



1430 West McCoy Lane
Santa Maria, CA 93455



1430 West McCoy Lane
Santa Maria, CA 93455