

# QUICKSLIDE™

Automated Slide Stainers





## Cat. no GPI

### ELECTRO-OPTICAL PROCESS

The GramPRO 1™ automatically adjusts for sample thickness, has individual preference controls for color intensity and is operator friendly.

### FAST, HANDS-FREE STAINING

Totally hands free, each slide is automatically stained, decolorized and counter stained in about 4 to 5 minutes. The GramPRO 1™ has an automatic cleaning cycle, has internal diagnostic software.

### SMALL COMPACT SIZE

The GramPRO 1™ is small enough to fit on any lab counter top.

### MINIMAL DAILY MAINTENANCE

Easy and fast. It's as simple as the push of a button.

## WHAT ARE THE BENEFITS OF THE GramPRO 1™ UNIT OVER MANUAL GRAM STAINING?

- The GramPRO 1™ is not "technique dependent". You get consistency no matter who processes the slide. With our patented electronic eye, decolorization is perfect every time.



## FEATURES

- Patented Decolorization process that optically monitors run-off fluid at a rate of 20 times per second
- Totally hands free process – Soft touch keypad
- Stain times in about 4 to 5 minutes
- Adjustable stain settings for color intensity
- Automatic calibration of each sample
- Adjustable sleep mode for energy efficiency
- Stain resistant metal cabinet fits on counter top
- Consistency for all sample types
- Automatic cleaning cycle
- Internal diagnostic software
- Audible on/off alarm to signal end of test
- Decolorize only cycle
- Manufactured in USA

- 
- There is no cross-contamination. The GramPRO 1™ uses fresh stain for each sample slide.
  - No more filtering of stains as with the manual process. This saves time and eliminates the mess.
  - The GramPRO 1™ has independently adjustable stain times on both the long and short cycles. It also has the option of decolorizing only if the operator chooses to do so.
  - The GramPRO 1™ allows the user to do other tasks while the slide is processing.
  - No technical expertise is required to do a quality Gram Stain. As long as the sample is properly heat fixed, anyone can perform the Gram Staining process with confidence.
  - Fewer fluids are used saving the lab disposal fees. The GramPRO 1™ does not have to be placed by sink – it can be placed beside a microscope.
  - Spend more time testing and less time cleaning your machines. Check out the competitors' recommended scheduled maintenance and you will appreciate the savings being offered by the GramPRO 1™.

**TOP SELLER!**



## Cat. no QS-001

The HemaPRO™ Slide Stainer represents a significant advancement in automated slide staining and is a valuable addition to any laboratory, whether it is the primary stainer used by small hospitals, physicians labs, stat labs, or used as a stat/backup instrument in larger labs. The HemaPRO™ can also be utilized in veterinarian labs for a variety of biological specimens.

### NEW FEATURE!

In addition to the soft agitation before the Rinse is drained, HemaPRO™ has an added soft agitation feature that occurs before the Buffer is drained. Through extensive testing we have determined this added feature puts off a cleaner slide with the same great staining characteristics of our modified Wright Giemsa stain formula. The unit needs no daily maintenance and only requires a tubing change once or twice annually.



## FEATURES

- State-of-the-art microprocessor, increased speed and ease of operator use
- Operator prompt LCD display
- Soft touch color-coded keypad
- Adjustable stain and buffer timing – so important for consistent lab results with multiple users
- Adjustable fill level
- “Cycles Remaining” display – computer monitors reagent and prompts when to reorder
- Injection molded case with built-in reagent platform
- Stainless steel reagent cannulas – oxidation resistant, increase system durability
- Color-coded reagent tubing for ease of operation
- Special formulated reagent kit – 125 tests per kit in one self-contained package

### SMALL COMPACT SIZE

The HemaPRO™ is one of the smallest automated slide stainers available. It requires only one square foot of counter space. In addition to its small size, the HemaPRO™ reagent system is included in the one square foot.

### FAST STAINING

The HemaPRO™ can stain a blood smear in as little as 60 seconds. In addition to the rapid staining, the HemaPRO™ eliminates cross contamination seen in the dip method and yields consistent results no matter who stains the slide.

### PERFORMS MULTIPLE STAINING PROCEDURES

In addition to staining peripheral blood smears, the HemaPRO™ has a special program that allows it to stain bone marrow samples. The HemaPRO also allows operator to adjust both stain and buffer times.



QUICKSLIDE™

1430 West McCoy Lane  
Santa Maria, CA 93455