

READYRACK™

A convenient way to bring media into a controlled area while avoiding outside contamination.

80 tube ReadyRack™

Convenient hand-holds



Remove ReadyRack™ from box outside of cleanroom.



Take into cleanroom, then remove plastic bag.



Spray down entire ReadyRack™ with isopropyl alcohol.

Packaged in a Clean Environment to Reduce Particulates

Tryptic Soy Broth (TSB), USP

16x125mm tube, 10ml fill, 80/pk..... K82BX

Fluid Thioglycollate Medium (FTM), USP

16x125mm tube, 10ml fill, 80/pk..... K212BX

Other tube formats and fill volumes are available upon request. Call (800) 266-2222 ext. 5696



Online Ordering Available

www.HardyDiagnostics.com/ReadyRack.html

HARDY
DIAGNOSTICS

READYRACK™



HARDY
DIAGNOSTICS

1430 West McCoy Lane
Santa Maria, CA 93455

I would like more information on this product.

Please fill in your information below and fax
the completed card to us at 805.346.2760.

Name _____

Title _____

Company _____

Address _____

City _____

State/Zip _____

Phone _____

Email _____

For more information on Ready Rack™ or any other Hardy Diagnostic Product:

Call toll free at 800.266.2222

Email us at Sales@HardyDiagnostics.com

Visit us online at www.HardyDiagnostics.com

A Culture of Service™