




TransPRO<sup>™</sup> TRANSPORT SYSTEM<sup>™</sup> CVM



**HARDY**  
DIAGNOSTICS  
A Culture of Service<sup>™</sup>



Online Ordering Available



102313TR

# Chlamydia, Viral, Mycoplasma Collection & Transport

TransPRO™ CVM Transport System is a specimen collection kit featuring PurFlock® Ultra swab and room temperature stability.

- Ready-to-use transport system
- Advanced kit design for reliability and ease-of-use
- Soft tip for patient comfort
- Flexible shaft
- Premium medical grade plastic components to improve product shelf life
- Stable medium at room temperature (see below)
- Economically priced

TransPRO™ CVM Transport System is a specimen collection kit featuring the PurFlock Ultra® Flocked swab and transport system. It is the ideal choice with excellent absorption and elution of specimens. For collection and preservation of Virus, Chlamydia, Mycoplasma, and Ureaplasma.

## Ordering Information

### Specimen Collection Kit with PurFlock® Ultra Swab

Adult/Pediatric,  
3ml, 50/pk.....TPCVM1




## A Comparison of Features

### Main Competitive Advantages in Favor of TransPRO™ CVM

TransPRO™ CVM		UTM™
Pack Size	50/pk	300/pk
Intended Use	Viruses, chlamydia, mycoplasma, ureaplasma	Viruses, chlamydia, mycoplasma, ureaplasma
Storage Temperature	2-25 °C	2-25 °C
Transport Temperature	4-25 °C for 48 hours	4-25 °C for 48 hours
Swab Type:	Naso-pharyngeal	Naso-pharyngeal
Swab-Cap Capture	Yes	Standard tip only
Swab Tip	PurFlock® Ultra (Patented New Technology)	Nylon Flocked (Patented)
Swab Tip Color	White	Yellow tinge
Cap Color	Red	Red
Contains Gelatin	Yes	Yes
Manufacture Location	USA	Italy
Shelf Life	18 months	12 months
Fill Volume	3ml	3ml
Cat. no.	TPCVM1	346C

\* UTM is a Trademark of Copan



**Swab-Cap Capture Feature:**

After the sample is collected, the swab is placed into the CVM Medium and the shaft is snapped at the score point. The lid is then screwed on. When the lid is unscrewed, the swab is connected to the lid.

Flexible Polystyrene Shaft



PurFlock® Ultra Flocked Swab Tip Technology for maximum specimen recovery. Clean, white, effective!

Vial:  
Overall Length  
16 x 100 mm  
(Compatible with  
most automation)

TransPRO™ CVM Medium



**Packaging includes:**

Vial with 3ml CVM Medium  
and one PurFlock® Ultra Flocked  
Swab, individually wrapped.

**Call 800.266.2222 for a  
free evaluation sample!**



**Made in the USA**





*A Culture of Service™*

1430 West McCoy Lane

Santa Maria, CA 93455

[www.HardyDiagnostics.com](http://www.HardyDiagnostics.com) | [Sales@HardyDiagnostics.com](mailto:Sales@HardyDiagnostics.com)

Distribution Centers: California | Washington | Utah | Arizona

Ohio | Texas | New York | Florida | North Carolina