



Water Testing for E. coli and Coliforms

colitag™

 **HARDY**
DIAGNOSTICS
A Culture of Service™



Online Ordering Available

02313TR

Water Testing for *E. coli* and Coliforms

What makes Colitag™ the most reliable and easy-to-use test?

- Fast results in as little as 16 hours
- Enhanced blue fluorescence to eliminate questionable positive *E. coli* results
- US EPA approved for MPN
- Detects one CFU of total coliform or *E. coli* bacteria per 100ml water sample
- Resuscitates chlorine-injured bacteria
- Growth inhibitors to provide selectivity without affecting recovery of target organisms

Order


| | |
|--|-----------|
| Colitag™ P/A - powdered media packet, 100/pk | 46000013 |
| Colitag™ Sample Container - with quick-dissolve sodium thiosulfate, screw cap with tamper evident seal, 100/pk | 46000121C |
| Colitag™ Comparator - Positive control for visual and fluorescent read-out, Each | 46000040 |
| UV Lamp - 6 watt, 365nm, long wave, hand held, Each | UVL56 |
| Colitag™ QC Standard Set - LyfoDisk™ Vial, 6 pellets per vial (2 vials/kit), Each | 5181L |
| Contains: <ul style="list-style-type: none">• <i>Escherichia coli</i> ATCC® 25922™• <i>Klebsiella pneumoniae</i> ATCC® 13883™• <i>Salmonella enterica</i> ATCC® 14028™ | |



Now approved for 16 hour read-out!



Add powder to water sample



Incubate for 16 hours



Positive for Coliform (yellow)



Positive for *E. coli* under UV light (fluorescent blue)



Call or scan QR code for a free sample and price quote! 800.266.2222



A Culture of Service™

1430 West McCoy Lane

Santa Maria, CA 93455

www.HardyDiagnostics.com | Sales@HardyDiagnostics.com

Distribution Centers: California | Washington | Utah | Arizona
Ohio | Texas | New York | Florida | North Carolina