

# Accuracy & Traceability have never been easier or more affordable.



**Introducing EZ-Accu Shot™, a quantitative microorganism preparation for Growth Promotion Testing.**

USP recommended strains • Self-contained kit • Precise CFU count • Delivers less than 100 CFU per 0.1 mL  
• 8-hour stability • Easy storage • No dry ice • Reduced shipping cost • 1 pellet per vial • No forceps required  
• Instant dissolve • Little to no prep time • No pre-incubation • Peel-off Certificate of Assay • On-line Certificate of Analysis • Cost-effective • Authentic • Traceable



phone: (800) 266-2222  
fax: (805) 614-9274  
website: [www.HardyDiagnostics.com](http://www.HardyDiagnostics.com)  
email: [Sales@HardyDiagnostics.com](mailto:Sales@HardyDiagnostics.com)

Santa Maria, CA  
Salt Lake City, UT  
Phoenix, AZ  
Centralia, WA  
Lima, OH  
Dallas, TX



*A culture of service ...*



## EZ- Accu Shot™ Kit Includes:

- 5 vials of a single quantitative microorganism (1 lyophilized pellet per vial)
- 5 vials of hydrating fluid (1.2 mL in each vial)
- Peel-off Certificate of Assay



Features	Benefits
Kit contains microorganism and hydrating fluid	Simplifies the ordering process, saving time and money
Instant dissolve, no pre-incubation	Little to no preparation time
One pellet per via	Eliminates the need for forceps, reducing the risk of contamination
Lower price than competitive products	Cost effective & economical
Eight hour stability	Extended time & flexibility to perform tests
Precision in CFU count	Increased customer confidence
Does not require freezing storage condition	Easy storage & no dry ice required, reducing shipping costs
Peel-off Certificate of Assay Label	Complete & convenient documentation
Online Certificate of Analysis	Documentation of phenotypic, macroscopic and microscopic test results
Authentic ATCC Licensed Derivative® Strains	USP recommended & traceable for assured quality
Product Warranty	Increased customer confidence

### Available Microorganism Strains:

Catalog #	Description
0317A	Clostridium sporogenes ATCC® 19404™**
0392A	Aspergillus brasiliensis ATCC® 16404™**
0443A	Candida albicans ATCC® 10231™**
0483A	Escherichia coli ATCC® 8739™**
0484A	Pseudomonas aeruginosa ATCC® 9027™**
0485A	Staphylococcus aureus subsp. aureus ATCC® 6538™**
0486A	Bacillus subtilis subsp. spizizenii ATCC® 6633™**
0487A	Clostridium sporogenes ATCC® 11437™**



phone: 800.266.2222  
 fax: 805.614.9274  
 www.HardyDiagnostics.com  
 Sales@HardyDiagnostics.com

Santa Maria, CA  
 Salt Lake City, UT  
 Phoenix, AZ  
 Centralia, WA  
 Lima, OH  
 Dallas, TX

EZ Accushot™ is manufactured by Microbiologics.

\* The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. MicroBiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

