

EnviroTrans™

*A ready-to-use swab & tube
pre-filled with your choice
of five solutions to fit
your specific needs!*



Collect and transport your environmental sample with efficiency and ease. Simply remove the cap, which is attached to the swab. Collect your sample by swabbing the surface, recap, and send to the lab.

A Self Contained System That is Easy to Use

The Hardy Diagnostics EnviroTrans™ is a ready-to-use 15 x 75mm polypropylene tube pre-filled with 5mls* of either Neutralizing Buffer, Lethen Broth, 0.85% Saline, DE Neutralizing Broth, or Neutralizing Saline. A dacron swab is affixed to a screw cap for easy handling. The entire unit is terminally sterilized. Excellent for detection or enumeration of microorganisms from environmental surfaces or equipment. Ideal for your HACCP program!

**The SRK05 is a .85% Saline solution with a .5ml fill.*

Call our Technical Service department for a recommendation to fit your specific needs.

How to Use the EnviroTrans™

When sampling a surface, the cap is removed and the moistened swab is rubbed across the sampling area. The cap is replaced on to the tube and the swab is transported to the laboratory for analysis.

Solutions to Fill Your Needs

Depending on your laboratory protocol you may choose from the following five solutions:

- ◆ **Neutralizing Buffer** for the detection of microorganisms found on surfaces disinfected with chlorine or quaternary ammonia compounds.
- ◆ **Lethen Broth** contains casein peptone and beef extract which promotes growth. Lecithin and Tween® are effective in neutralizing phenols, hexachlorophene, formalin, and ethanol. Great for use with Petrifilm™.
- ◆ **0.85% Saline** is a sterile osmotically neutral transport solution.
- ◆ **DE Neutralizing Broth** used to neutralize antiseptics and disinfectants as well as detect organisms remaining after treatment.
- ◆ **Neutralizing Saline** neutralizes the broadest spectrum of antiseptics and disinfectant chemicals while facilitating survival of the organism in the sample during transport. Recommended for enumeration.

EnviroTrans™ with Neutralizing Buffer	20/box	SRK15
EnviroTrans™ with Lethen Broth	20/box	SRK25
EnviroTrans™ with 0.85% Saline	20/box	SRK35
EnviroTrans™ with 0.85% Saline, .5ml	20/box	SRK05
EnviroTrans™ with DE Neutralizing Broth	20/box	SRK45
EnviroTrans™ with Neutralizing Saline	20/box	SRK55



Santa Maria, CA 93455
Salt Lake City, UT 84070
Phoenix, AZ 85210

Tel: (800) 266-2222
Fax: (805) 346-2760

Website: hardydiagnostics.com
Email: sales@hardydiagnostics.com