

StrepPRO™ Grouping Kit

A unique streptococcal grouping latex agglutination kit for the rapid identification of group A, B, C, D, F, and G *Streptococcus* from culture.



StrepPRO™ Grouping Kit

60 Test Kit..... PL030HD

Kit components are also available separately.



Positive Reaction



Negative Reaction

StrepPRO™ Grouping Kit

- Group A.... PL031HD
- Group B.... PL032HD
- Group C..... PL033HD
- Group G..... PL036HD

- Extraction Reagent 1..... PL031HD
- Extraction Reagent 2..... PL033HD
- Extraction Reagent 3..... PL034HD

- Positive Control..... PL035HD
- Testing Cards..... PL036HD

Faster

Accurate results in one minute at room temperature. Utilizes modified acid extraction technology. No incubation required.

More Convenient

All color-coded reagents are ready-to-use and remain in one space-saving workstation. 60 tests/kit
*Kit includes Group D Latex.

Easier to Read

Blue latex particles on a white reaction card with eight reaction circles creates an exceptionally easy-to-read result.

Greater Sensitivity

Greater than 98% sensitive with as little as one colony required.

Compact size conserves valuable refrigerator space

 **HARDY**
DIAGNOSTICS
A Culture of Service™

StreppRO™ Grouping Kit™

I would like more information on these products.
Please fill in your information below and mail it or snap a photo of this postcard and send it to Sales@HardyDiagnostics.com

Name _____
Title _____
Company _____
Address _____
City _____
State/Zip _____
Phone _____
Email _____

**For more information and ordering
StreppRO™ and other
Hardy products visit our website at
www.HardyDiagnostics.com
or scan our QR code.**

Call toll free at 800.266.2222 ext. 5696
Email us at Sales@HardyDiagnostics.com



HARDY
DIAGNOSTICS
A Culture of Service™
1430 West McCoy Lane
Santa Maria, CA 93455