

Water Sampling Vials

*Collect your
water sample and
neutralize all chlorine!*



A water sampling vial that meets all EPA requirements for the microbiological water and wastewater testing!

Meets all EPA Requirements.

A 120ml, pre-sterilized, hinged cap vial with or without a sodium thiosulfate tablet. Meets all EPA requirements for a microbiological container for water or wastewater testing. Compatible with Idexx's Colilert™ method of testing for coliforms and *E. coli*.

Important Features!

- ◆ 120ml total vial capacity.
- ◆ 100ml raised sample fill ring.
- ◆ Durable, contact clarity polypropylene.
- ◆ Pre-sterilized with or without Sodium Thiosulfate tablet.
- ◆ Innovative, snap-down, Lock-N-Seal latch, eliminates accidental opening during transport. Custody tie, ensures that the samples have not been tampered with.

EPA Vial	With Tablet	200/case	EPA12083P
EPA Vial	Without Tablet	200/case	EPA12083

HIDx



Santa Maria, CA 93455
Salt Lake City, UT 84070
Phoenix, AZ 85210

Tel: (800) 266-2222
Fax: (805) 346-2760

Website: hardydiagnostics.com
Email: sales@hardydiagnostics.com